

By **YORCO**
Scientific Udyog Pvt. Ltd.

EVERY PACKAGE
BEING STERILIZED
SAVES LIFE



PRODUCT PORTFOLIO

ABOUT US

STERI has been at the forefront of manufacturing sterilization and infection control products for over 30 years. Being a family owned business, we have developed a reputation for our quality and reliability. At **STERI**, we believe in providing excellent products and exceptional services to our clients. Our experience with **over 50,000** installations and strive for excellence has resulted in **STERI** being a pioneer in the field.

At **STERI** we ensure exceptional customer support by providing expert **pre and post sales support services**, to satisfy customer expectations. Our installation team provides in depth one on **one training and demos and our sales and technical team is highly skilled** to trouble shoot any technical issues online and in person nationwide. We highly value customer feedback which contributes to the continuous improvement of our products and support services.

WHY US

CERTIFICATIONS



ISI

Indian Standard Institute



FDA

The Food and Drug Administration



CE - EU

Standards for health, safety and environmental protection



GMP

Good Manufacturing Practice



ISO:9001

Quality Managements System



ISO 13485

Medical device quality management system



IEC 61010-2

Safety requirements for electrical equipment for measurement, control, and laboratory use



IEC 60601-1-2016

Technical standard for the safety and essential performance of medical electrical equipment

SERVICES




Pre - Sale Service



Post - Sale Service



1 Year Warranty



Free Onsite Training



Fast Shipping



Free Installation

QUALITY



Hydraulically Tested at 2.5x Pressure



Multiple Quality Check

Sterilizer Series

STEAM STERILIZER



YSU-404 STERI VERTICAL STEAM STERILIZER

- Automatic Microprocessor Touch Screen, Semi-Automatic and Manual Control Systems

| | |
|--------------|--------------|
| Capacity | 35 L, 75 L |
| Working Temp | 121°C, 134°C |



YSU-405 STERI HORIZONTAL STEAM STERILIZER

- Semi-Automatic and Manual Control Systems

| | |
|--------------|---|
| Capacity | 75 L, 138 L, 176 L, 235 L, 254 L, 485 L |
| Working Temp | 121°C, 134°C |



YSU-608 STERI AUTOMATIC STEAM STERILIZER

- Available in single, double and sliding door
- Microprocessor Touch Screen Control System

| | |
|--------------|---|
| Capacity | 75 L, 138 L, 176 L, 235 L, 254 L, 485 L |
| Working Temp | 100°C-134°C |



YSU-617 TABLE TOP PORTABLE AUTOCLAVE

- Automatic Microprocessor Touch Screen and Manual Control Systems with sliding doors

| | |
|--------------|-------------|
| Capacity | 14L, 19L |
| Working Temp | 100°C-134°C |

GAS STERILIZER



YSU-606 STERI ETO HORIZONTAL STERILIZER

- Automatic Microprocessor Touch Screen Control System

| | |
|------------------------|------------------------|
| Capacity (Rectangular) | 176L, 235L, 254L, 485L |
| Capacity (Cylindrical) | 75 L, 138 L |
| Working Temp | 37°C- 60°C |



YSU-406 STERI ETO TABLE TOP STERILIZER

- Automatic Microprocessor Touch Screen Control System

| | |
|------------------------|--------------|
| Capacity (Rectangular) | 176 L, 235 L |
| Capacity (Cylindrical) | 35 L, 75 L |
| Working Temp | 37°C- 60°C |



YSU-607 STERI FORMALDEHYDE STERILIZER

- Automatic Microprocessor Touch Screen Control System

| | |
|------------------------|------------------------|
| Capacity (Rectangular) | 176L, 235L, 254L, 485L |
| Capacity (Cylindrical) | 75 L, 138 L |
| Working Temp | 60°C- 80°C |

Effective against vegetative bacteria, Mycobacteria, B. Atrophaeus and G. Stearothermophilus Spores and Candida Albicans

OTHER STERILIZERS



YSU-701 STERI UV STERILIZER



YSU-623 STERI DRY STERILIZER



YSU-409 STERI HOT & COLD WATER STERILIZER



YSU-613 STERI BOWL STERILIZER



YSU-900 STERI INSTRUMENT STERILIZER

Decontamination Series



YSU-609 STERI DISINFECTOMAT WASHER DISINFECTOR

- Automatic Microprocessor Touch Screen Control System
- 6 washing programs, auto detergent water heating and water discharging system

YSU-615 STERI SCRUB STATION

- Proximity sensor for Automatic Operation
- Available as wall mounted or freestanding units
- High grade 304 construction
- 1,2 and 3 sink configuration

YSU-1150 STERI HAND WASH SINK

- High grade SS-304 construction
- 1,2 and 3 sink configuration

Distillation Series

YSU-410 STERI SINGLE DISTILLATION WATER STILL



- Complete SS-304 construction
- Produces pyrogen free distilled water

| | |
|-----------------|--|
| | Distillate L/hr |
| Capacity | 5 L/hr, 10L/hr, 20 L/hr, 40L/hr, 60 L/hr, 100 L/hr |

YSU-410A STERI DOUBLE DISTILLATION WATER STILL



- Complete SS-304 construction
- Produces pyrogen free distilled water

| | |
|-----------------|---------------------------------|
| | Distillate L/hr |
| Capacity | 5 L/hr, 10L/hr, 20 L/hr, 40L/hr |

Bioengineering Series

YSU-551 STERI FERMENTER MANUAL



- Manual and Automatic Advanced Microprocessor PLC Touch Screen Control System
- Also available with a Inbuilt sterilization

| | |
|---------------------|----------------------|
| Capacity | 10L, 50L, 100L, 200L |
| Working Temp | 37°C- 60°C |

YSU-555 STERI AUTOMATIC BABY FERMENTER



- All functions are controlled via Automatic Advanced Microprocessor PLC Touch Screen Control System
- Also available with a Glass vessel

| | |
|---------------------|------------|
| Capacity | 5L |
| Working Temp | 37°C- 60°C |

Surgical Series

YSU-1129 STERI DIGITAL LED LIGHT SOURCE



STERI in collaboration with ARK has developed a series of high feature low cost medical led light source for Laparoscopy , Endoscopy , ENT , Dental etc.

- 100 Watt LED at 6000° – 6500°K ensures constant daylight colour rendition
- Novel heat protection filter, analog and digital electronic light intensity regulator

Waste Management Series

YSU-611 CYLINDRICAL MEDICAL WASTE STERILIZER WITH SHREDDER



- Advanced Microprocessor PLC Touch Screen Control System
- Also available with attached shredder

| | |
|-------------------------------|------------------------------------|
| Capacity (Rectangular) | 432 L, 648 L, 810 L, 972 L, 1215 L |
| Capacity (Cylindrical) | 235 L, 340 L, 530 L |

YSU-612 WASTE SHREDDER



- Designed to shred plastic waste such as blood bags, transfusion sets, plastic syringes, scalpel, blades etc

| | |
|-----------------|-----------------------------------|
| Capacity | 10 kg, 20 kg, 30 kg, 40 kg, 50 kg |
|-----------------|-----------------------------------|



LINEN TROLLEY



WASTE SEGREGATION TROLLEY



BIO HAZARD TROLLEY

Furniture & CSSD Items

1. Packing Table
2. Paper Table
3. Dry Work Table
4. SS Table
5. Gauze Cutting Table
6. CSSD Rack
7. Basket Rack



1. Semi Fowler Bed
2. Fowler Bed
3. Hospital bed
4. ICU Bed
5. Bed Side Locker
6. Examination Couch
7. SS Almera
8. Baby Bassinet



Our Esteemed Clients

Get In Touch

HOSPITALS AND MEDICAL COLLEGES



GOVERNMENT INSTITUTION



PHARMACEUTICAL LAB



FOOD AND BEVERAGE



AGRICULTURAL LAB



MARKETED BY



📍 :D-149 Preet Vihar,
New Delhi 110092.

✉ : ysuptld@gmail.com

🌐 : yorcoscientific.com

☎ : 8595293382

MANUFACTURED BY



YORCO
Scientific Udyog Pvt.Ltd

📍 : 13 Loni road, site 2,
Mohan nagar,
Ghaziabad,
Uttar Pradesh 201007.

✉ : ysuptld@gmail.com

🌐 : yorcoscientific.com

☎ : 9350000607