



## YSU-409 Steri Hot & Cold Water

**Sterilizer as per IS: 7455**

Usage/Application	Hospital
Type Of Sterilizer	Hot And Cold Water Sterilizer
Material	Stainless Steel
Insulation Wall	Triple Wall
Automation Grade	Fully Automatic
Brand	Steri
Size	30 to 50 Liters
Power Supply	12 Kw



### **-: DESCRIPTION :-**

The unit is mounted on a sturdy M.S. stand, consisting of two separate st. steel tank one is for **sterile hot water** and other is for **sterile cold water**.

The unit comprises of individual control such as **safety valve, pressure gauge, steam trap, self sterilizing water filter, water level indicator, air filter and draw off water valve with socket.**

Each tank has an individual heater bank, **low water cut off, pressure controller & contactor.**  
The unit is made as per IS standards i.e. IS: 7455 (1974).

**\*Useful in operation theaters where quick sterile hot or cold water is required.**

**Yorco Scientific Udyog Pvt. Ltd.**

Add-13, Loni Road, Site-2, Mohan Nagar, Ghaziabad (UP)-201007.

ysupltd@gmail.com / www.yorcoscientific.com / Ph.-91-935000607, 8595293382.