

By **YORCO**
Scientific Udyog Pvt. Ltd.

EVERY PACKAGE
BEING STERILIZED
SAVES LIFE



PRODUCT PORTFOLIO

Message from CEO



Mr Suresh Arora

Thank you for taking out the time to go through our product brochure. My name is Suresh Arora, the Director and CEO at Yorco Scientific Udyog Private limited. We have over 60 years of experience of manufacturing medical devices and completing over 50000 installation throughout India and abroad. Through these 60 years, you have immensely helped us grow through your valuable feedback. Learning from your inputs we are able to provide finest quality equipment's with upgraded features, improvised solutions and prompt after sale service.

While pursuing our vision to serve humanity by providing high-quality medical equipment for global wellbeing, we are fully committed to strengthen our capabilities and to further evolve the organization through your invaluable feedback. We have undergone a remarkable evolution with initiatives such as Innovative Product Development, Introducing automation like CNC machining, Lean Manufacturing, Online servicing, streamlining processes. These initiatives have helped us to become one of the leading Indian Manufacturers for scientific equipment. Thank you for letting us get better in order to serve you better.

Our Team



Kush Arora
Chief Operating Officer
Biomedical Engineer (M.S. B.Tech)



Kuldeep Kumar Agarwal
Production Manager
Mechanical Engineer (B.Tech)



Pankaj Kumar
Quality Manager
Mechanical Engineer (B.Tech)



ABOUT US

STERI has been at the forefront of manufacturing Sterilization and Infection control products for over **60 years**. With wide business base, complete in-house design and manufacturing facility, we have developed a reputation for our quality and reliability. At Steri we provide excellent products and exceptional customer service to our clients. Our experience with **over 50,000** installations and strive for excellence has resulted in **STERI** being a pioneer in the field.

At **STERI** we ensure exceptional customer support by providing expert **pre-sale advise, technical consultation, site inspection and after sales support services** to satisfy customer. Our installation team provides in depth **one on one training and demos** during installation. Our **sales and technical team is highly skilled** to trouble shoot any rare technical issues online and in person nationwide. We highly value customer feedback which contributes to the continuous improvement of our products and support services.

WHY US

CERTIFICATIONS



ISI

Indian Standard Institute



FDA

The Food and Drug Administration



CE - EU

Standards for health, safety and environmental protection



GMP

Good Manufacturing Practice



ISO:9001

Quality Managements System



ISO 13485

Medical device quality management system



IEC 61010-2

Safety requirements for electrical equipment for measurement, control, and laboratory use



IEC 60601-1-2016

Technical standard for the safety and essential performance of medical electrical equipment

SERVICES




Pre - Sale Service



After - Sale Service



1 Year Warranty



Free Onsite Training



Fast Shipping



Free Installation

QUALITY



Hydraulically Tested at 2.5x Pressure



Multiple Quality Check

Sterilizer Series

STEAM STERILIZER



YSU-404 STERI VERTICAL STEAM STERILIZER

- Automatic Microprocessor Touch Screen, Semi-Automatic and Manual Control Systems

Capacity	35 L, 75 L
Working Temp	121°C, 134°C



YSU-405 STERI HORIZONTAL STEAM STERILIZER

- Semi-Automatic and Manual Control Systems

Capacity	75 L, 138 L, 176 L, 235 L, 254 L, 485 L
Working Temp	121°C, 134°C



YSU-608 STERI AUTOMATIC STEAM STERILIZER

- Available in single, double and sliding door
- Microprocessor Touch Screen Control System

Capacity	75 L, 138 L, 176 L, 235 L, 254 L, 485 L
Working Temp	100°C-134°C



YSU-617 TABLE TOP PORTABLE AUTOCLAVE

- Automatic Microprocessor Touch Screen and Manual Control Systems with sliding doors

Capacity	14L, 19L
Working Temp	100°C-134°C

GAS STERILIZER



YSU-606 STERI ETO HORIZONTAL STERILIZER

- Automatic Microprocessor Touch Screen Control System

Capacity (Rectangular)	176L, 235L, 254L, 485L
Capacity (Cylindrical)	75 L, 138 L
Working Temp	37°C- 60°C



YSU-406 STERI ETO TABLE TOP STERILIZER

- Automatic Microprocessor Touch Screen Control System

Capacity (Rectangular)	176 L, 235 L
Capacity (Cylindrical)	35 L, 75 L
Working Temp	37°C- 60°C



YSU-607 STERI FORMALDEHYDE STERILIZER

- Automatic Microprocessor Touch Screen Control System

Capacity (Rectangular)	176L, 235L, 254L, 485L
Capacity (Cylindrical)	75 L, 138 L
Working Temp	60°C- 80°C

Effective against vegetative bacteria, Mycobacteria, B. Atrophaeus and G. Stearotherophilus Spores and Candida Albicans

OTHER STERILIZERS



YSU-701 STERI UV STERILIZER



YSU-623 STERI DRY STERILIZER



YSU-409 STERI HOT & COLD WATER STERILIZER



YSU-613 STERI BOWL STERILIZER



YSU-900 STERI INSTRUMENT STERILIZER

Decontamination Series



YSU-609 STERI DISINFECTOMAT WASHER DISINFECTOR

- Automatic Microprocessor Touch Screen Control System
- 6 washing programs, auto detergent water heating and water discharging system



YSU-615 STERI SCRUB STATION

- Proximity sensor for Automatic Operation
- Available as wall mounted or freestanding units
- High grade 304 construction
- 1, 2 and 3 sink configuration



YSU-1150 STERI HAND WASH SINK

- High grade SS-304 construction
- 1, 2 and 3 sink configuration

Distillation Series

YSU-410 STERI SINGLE DISTILLATION WATER STILL



- Complete SS-304 construction
- Produces pyrogen free distilled water

	Distillate L/hr
Capacity	5 L/hr, 10L/hr, 20 L/hr, 40L/hr, 60 L/hr, 100 L/hr

YSU-410A STERI DOUBLE DISTILLATION WATER STILL



- Complete SS-304 construction
- Produces pyrogen free distilled water

	Distillate L/hr
Capacity	5 L/hr, 10L/hr, 20 L/hr, 40L/hr,

Bioengineering Series

YSU-551 STERI FERMENTER MANUAL



- Manual and Automatic Advanced Microprocessor PLC Touch Screen Control System
- Also available with a Inbuilt sterilization

Capacity	10L, 50L, 100L, 200L
Working Temp	37°C- 60°C

YSU-555 STERI AUTOMATIC BABY FERMENTER



- All functions are controlled via Automatic Advanced Microprocessor PLC Touch Screen Control System
- Also available with a Glass vessel

Capacity	5L
Working Temp	37°C- 60°C

Surgical Series

YSU-1129 STERI DIGITAL LED LIGHT SOURCE



STERI in collaboration with ARK has developed a series of high feature low cost medical led light source for Laparoscopy , Endoscopy , ENT , Dental etc.

- 100 Watt LED at 6000° – 6500°K ensures constant daylight colour rendition
- Novel heat protection filter, analog and digital electronic light intensity regulator

Waste Management Series

YSU-611 CYLINDRICAL MEDICAL WASTE STERILIZER WITH SHREDDER



- Advanced Microprocessor PLC Touch Screen Control System
- Also available with attached shredder

Capacity (Rectangular)	432 L, 648 L, 810 L, 972 L, 1215 L
Capacity (Cylindrical)	235 L, 340 L, 530 L

YSU-612 WASTE SHREDDER



- Designed to shred plastic waste such as blood bags, transfusion sets, plastic syringes, scalpel, blades etc

Capacity	10 kg, 20 kg, 30 kg, 40 kg, 50 kg
-----------------	-----------------------------------



LINEN TROLLEY



WASTE SEGREGATION TROLLEY



BIO HAZARD TROLLEY

Furniture & CSSD Items

1. Packing Table
2. Paper Table
3. Dry Work Table
4. SS Table
5. Gauze Cutting Table
6. CSSD Rack
7. Basket Rack



1. Semi Fowler Bed
2. Fowler Bed
3. Hospital bed
4. ICU Bed
5. Bed Side Locker
6. Examination Couch
7. SS Almira
8. Baby Bassinet



Our Approach



1

Understanding Customer Needs



2

Customization as Required



3

Place the Order




6

Laboratory Testing



5

Material Procurement




4

Developing Design Files



PRODUCTION
PRE



7

Fabrication




8

Hydraulic Testing



9

Multiply Quality Testings




10

Final Assembly



PRODUCTION




13

Packaging



12

Finishing



11

Factory Acceptance Testing



14

Fast Shipping



15

Installation



16

On-Site Training

STERI continues to deliver superior-quality products that are manufactured in accordance with current Good Manufacturing Practices (cGMP).



17

Guaranteed Customer Satisfaction

Our Esteemed Clients

HOSPITALS AND MEDICAL COLLEGES



GOVERNMENT INSTITUTION



PHARMACEUTICAL LAB



FOOD AND BEVERAGE



AGRICULTURAL LAB



Get In Touch

MARKETED BY



Arktech Scientific Pvt.Ltd

📍 :D-149 Preet Vihar
New Delhi 110092.

✉ :arkscientific.sales@gmail.com

🌐 :yorcoscientific.com

☎ :8595293382

MANUFACTURED BY



YORCO
Scientific Udyog Pvt.Ltd

📍 : 13 Loni road, site 2,
Mohan nagar, Ghaziabad,
Uttar Pradesh 201007.

✉ : ysuptld@gmail.com

🌐 : yorcoscientific.com

☎ :9350000607