



By **YORCO**  
Scientific Udyog Pvt. Ltd.

# SINGLE DISTILLATION WATER STILL

**STERI** has been at the forefront of manufacturing Sterilization and Infection control products for over **60 years**. With wide business base, complete in-house design and manufacturing facility, we have developed a reputation for our quality and reliability. At Steri we provide excellent products and exceptional customer service to our clients. Our experience with **over 50,000** installations and strive for excellence has resulted in **STERI** being a pioneer in the field.

### Rust-Resistant Construction

High-Quality Stainless Steel is used to avoid corrosion and contamination

### Air Tight Sealing

Jointless steam resistant silicon gasket for effective sealing and safety

### Unique Condenser System

The condenser design ensures that the droplets of distilled water remain in contact with the cooling coil for the longest possible time achieving the distillate temperature of 60°C

### High Quality Manufacturing

Manufactured conforming to IS: No.3830 and European CE standards and equipped with ISI marked water immersion heating elements



### High Purity Water

Water Distillation equipment provides pyrogen free distillate

### Narrow and Compact

The slim design allows the machine to fit in a small space for maximum area utilization

### Advanced Sensor System

Steri Water Distillation comes equipped with state of the art smart electronic water sensor for accurate detection and heater protection by controlling the water

### Exquisite Finish

Equipment is mirror polished internally and matt polished externally to provide a smooth finish without any sharp edges and crevice

## Capacity Available

|   | VOLUME         | LOAD  |
|---|----------------|-------|
| 1 | 5 Liter/ hour  | 4 KW  |
| 2 | 10 Liter/ hour | 8 KW  |
| 3 | 20 Liter/ hour | 15 KW |
| 4 | 40Litre/ hour  | 30 KW |

Voltage 220V, 50 HZ  
Phase AC

## Accessories

1. INCLUDED WITH THE PRODUCT
  - Supplied with cord and plug to work on 220/230 volts AC
  - Sturdy Steel Stand (Powder Coated)
2. AVAILABLE FOR PURCHASE
  - Heater element
  - Joint less silicon gasket



## Other Features

All units are made as per IS : 3830 (ISI), CE, ISO: 9001, and ISO:13485 standards

Each unit is fitted with ISI marked immersion heating elements

The instrument is equipped with automatic low water protection system and an easy to use electrical control box

The stand is powder coated for rust free operation and lustrous appearance.

Fuse and current leakage breaker is provided for safety in case of voltage fluctuation or short circuit.

All the electrical and plumbing fixtures are assembled in our factory for quick and easy installation on site

Single Distillation Water Still comes with one year manufacturer's warranty

### Customization of Products :

- Double-walled SS construction with insulation available upon request
- Construction in Stainless Steel 316 available on request
- Fully Automatic control system
- In-built water softener system

## Why Us ?



Affordable



Low Maintenance



Customizable



Water Saving



Hydraulically Tested at 2.5x Pressure



Multiple Quality Check




1 Year Warranty



Pre - Sale Service



After- Sale Service



Free Onsite Training



Fast Shipping



Free Installation

## Our Esteemed Clients



## Get In Touch



**STERI** - BY YORCO Scientific Udyog Pvt.Ltd

📍 : 13 Loni road, site 2,  
Mohan nagar, Ghaziabad,  
Uttar Pradesh 201007.

✉ : [ysupltd@gmail.com](mailto:ysupltd@gmail.com)

🌐 : [yorcoscientific.com](http://yorcoscientific.com)

☎ : 9350000607