

By **YORCO**
Scientific Udyog Pvt. Ltd.



EVERY PACKAGE BEING STERILIZED
SAVES LIFE



VERTICAL STERILIZER



About Us

STERI has been at the forefront of manufacturing Sterilization and Infection control products for over **60 years**. With wide business base, complete in-house design and manufacturing facility, we have developed a reputation for our quality and reliability. At Steri we provide excellent products and exceptional customer service to our clients.

Our experience with **over 50,000** installations and strive for excellence has resulted in **STERI** being a pioneer in the field. We not only cater to the hospital sector but have a diverse clientele including medical colleges, dental clinics, laboratories, veterinary hospitals, pharmaceutical companies, food & beverage industries and agricultural institutes.

Experience

Service

At **STERI** we ensure exceptional customer support by providing expert **pre-sale advise, technical consultation, site inspection and after sales support services** to satisfy customer. Our installation team provides in depth **one on one training and demos** during installation. Our **sales and technical team is highly skilled** to trouble shoot any rare technical issues online and in person nationwide. We highly value customer feedback which contributes to the continuous improvement of our products and support services.

Message from CEO



Mr Suresh Arora

Thank you for taking out the time to go through our product brochure. My name is Suresh Arora, the Director and CEO at Yorco Scientific Udyog Private limited. We have over 60 years of experience of manufacturing medical devices and completing over 50000 installation throughout India and abroad. Through these 60 years, you have immensely helped us grow through your valuable feedback. Learning from your inputs we are able to provide finest quality equipment's with upgraded features, improvised solutions and prompt after sale service.

While pursuing our vision to serve humanity by providing high-quality medical equipment for global wellbeing, we are fully committed to strengthen our capabilities and to further evolve the organization through your invaluable feedback. We have undergone a remarkable evolution with initiatives such as Innovative Product Development, Introducing automation like CNC machining, Lean Manufacturing, Online servicing, streamlining processes. These initiatives have helped us to become one of the leading Indian Manufacturers for scientific equipment. Thank you for letting us get better in order to serve you better.

Our Team



Kush Arora
Chief Operating Officer
Biomedical Engineer (M.S. B.Tech)



Kuldeep Kumar Agarwal
Production Manager
Mechanical Engineer (B.Tech)



Pankaj Kumar
Quality Manager
Mechanical Engineer (B.Tech)



Why Us ?

CERTIFICATIONS



ISI

Indian Standard Institute



FDA

The Food and Drug Administration



CE - EU

Standards for health, safety, and environmental protection



GMP

Good Manufacturing Practice



ISO:9001

Quality Managements System



ISO 13485

Medical device quality management system



IEC 61010-2


Safety requirements for electrical equipment for measurement, control, and laboratory use



IEC 60601-1-2016

Technical standard for the safety and essential performance of medical electrical equipment

DESIGN FEATURES



Uniform Heat Distribution




Rust Proof



Low Maintenance



Affordable



Compact Design



Customizable



User Friendly



Water Saving

SERVICES




Pre - Sale Service



After - Sale Service



1 Year Warranty



Free Onsite Training



Fast Shipping



Free Installation

QUALITY




Hydraulically Tested at 2.5x Pressure



Multiple Quality Check

Product Features



Preset Sterilization Cycles

- Instrument Cycle
- Fabric Sterilization Cycle
- Liquid Sterilization Cycle
- Wrapped Sterilization Cycle
- Customer Variable Cycle
- Leak Test Cycle

*More cycles can be added as per customer requirements

Accessories Often Bought Together



Sterilization Tape



SS Mesh Basket



Stackable
SS Mesh Basket



Dressing Drum







Basket Stand

Sizes Available

DIAMETER	X	DEPTH	VOLUME	LOAD
300mm	X	500mm	35 liters	4 kw
400mm	X	600mm	75 liters	6 kw

*Being an OEM we can always modify the chamber sizes according to customer requirements.

Product Range

	BRONZE	SILVER	GOLD	PLATINUM
OUR PRODUCTS	 YSU-404 B	 YSU-404 S	 YSU-404 G	 YSU-404 P
Body				
Chamber	SS304*	SS304*	SS304*	SS316*
Steam Jacket	SS304*	SS304*	SS304*	SS316*
Chamber Cover	SS304*	SS304*	SS304*	SS304*
Lid	SS304*	SS304*	SS304*	SS316*
Lid Locking	Wingnut	Wingnut	Wingnut	Radial
Outer Cover	●	●	●	✓ (MS/SS)
Other Features				
Piping SS304	✓	✓	✓	✓
Safety Valve & Vacuum Breaker	✓	✓	✓	✓
Low Water Cut Off	✓	✓	✓	✓
Moisture Trap	✓	✓	✓	✓
Dial Thermometer	✓	✓	✓	✓
Analog Pressure Gauge	✓	✓	✓	✓
Analog Compound Gauge	✓	✓	✓	✓
Automatic Pressure Control	✓	✓	✓	✓
Automation Device	●	DTC	PID	PLC
Automatic Ster Controlled Valves	●	✓	✓	✓
Digital Temperature Controller	●	✓	✓	✓
Digital Temperature Display	●	✓	✓	✓
Digital Timer	●	✓	✓	✓
Alarm	●	✓	✓	✓
Digital Pressure Display	●	●	✓	✓
Printer	●	●	optional	✓
Touch Screen Display	●	●	●	✓
Pre Set Sterilization Cycle	●	●	●	✓
Automatic Water feeding	●	●	●	✓
Auto Exhaust	●	●	●	✓

Customization of products can be done as per client's requirements

*Highly corrosion resistant SS316/ SS316L/ SS316TI are also available on request

Our Approach



1 Understanding Customer Needs



2 Customization as Required



3 Place the Order


PRODUCTION
PRE




4 Laboratory Testing



5 Material Procurement



6 Developing Design Files




7 Fabrication



8 Hydraulic Testing




9 Multiply Quality Testings



10 Final Assembly


PRODUCTION



11 Packaging



12 Finishing



13 Factory Acceptance Testing

POST
PRODUCTION



14 Fast Shipping



15 Installation



16 On-Site Training

STERI continues to deliver superior-quality products that are manufactured in accordance with current Good Manufacturing Practices (cGMP).



17 Guaranteed Customer Satisfaction

Our Esteemed Clients

Hospitals and Medical Colleges



Pharmaceutical Lab



Food and Beverage



Agricultural Labs



Get in Touch

Marketed by



Arktech Scientific Pvt.Ltd
:D-149 Preet Vihar
New Delhi 110092.
✉ :arkscientific.sales@gmail.com
🌐 :yorcoscientific.com
☎ :8595293382

Manufactured by



YORCO Scientific Udyog Pvt.Ltd
: 13 Loni road, site 2,
Mohan nagar, Ghaziabad,
Uttar Pradesh 201007.
✉ : ysuptld@gmail.com
🌐 : yorcoscientific.com
☎ :935000607