

By **YORCO**
Scientific Udyog Pvt. Ltd.



EVERY PACKAGE BEING STERILIZED
SAVES LIFE



VERTICAL STERILIZER



About Us

STERI has been at the forefront of manufacturing Sterilization and Infection control products for over **60 years**. With wide business base, complete in-house design and manufacturing facility, we have developed a reputation for our quality and reliability. At Steri we provide excellent products and exceptional customer service to our clients.

Our experience with **over 50,000** installations and strive for excellence has resulted in **STERI** being a pioneer in the field. We not only cater to the hospital sector but have a diverse clientele including medical colleges, dental clinics, laboratories, veterinary hospitals, pharmaceutical companies, food & beverage industries and agricultural institutes.

Experience

Service

At **STERI** we ensure exceptional customer support by providing expert **pre-sale advise, technical consultation, site inspection and after sales support services** to satisfy customer. Our installation team provides in depth **one on one training and demos** during installation. Our **sales and technical team is highly skilled** to trouble shoot any rare technical issues online and in person nationwide. We highly value customer feedback which contributes to the continuous improvement of our products and support services.

Why Us ?

CERTIFICATIONS



ISI
Indian Standard Institute



FDA
The Food and Drug Administration



CE - EU
Standards for health, safety, and environmental protection



GMP
Good Manufacturing Practice



ISO:9001
Quality Managements System



ISO 13485
Medical device quality management system




IEC 61010-2
Safety requirements for electrical equipment for measurement, control, and laboratory use



IEC 60601-1-2016
Technical standard for the safety and essential performance of medical electrical equipment

DESIGN FEATURES



Uniform Heat Distribution




Rust Proof




Low Maintenance




Afordable



Compact Design



Customizable



User Friendly



Water Saving

SERVICES




Pre - Sale Service



Post - Sale Service



1 Year Warranty



Free Onsite Training



Fast Shipping



Free Installation

QUALITY



Hydraulically Tested at 2.5x Pressure



Multiple Quality Check

Product Features



Preset Sterilization Cycles

- Instrument Cycle
- Linen Sterilization Cycle
- Liquid Sterilization Cycle
- Customer Variable Cycle
- Bowie Dick Test Cycle
- Leak Test Cycle

Sizes Available

DIAMETER	X	DEPTH	VOLUME	LOAD
300mm	X	500mm	35 liters	4 kw
400mm	X	600mm	75 liters	6 kw

Accessories Often Bought Together



Sterilization Tape



SS Mesh Basket



Stackable
SS Mesh Basket







Dressing Drum



Basket Stand

Product Range

	BRONZE	SILVER	GOLD	PLATINUM
OUR PRODUCTS	 YSU-404 B	 YSU-404 S	 YSU-404 G	 YSU-404 P
Body				
Chamber	SS304	SS304	SS304	SS316
Steam Jacket	SS304	SS304	SS304	SS316
Chamber Cover	SS304	SS304	SS304	SS304
Lid	SS304	SS304	SS304	SS316
Lid Locking	Wingnut	Wingnut	Wingnut	Radial
Outer Cover		●	●	✓ (MS/SS)
Other Features				
Piping SS304	✓	✓	✓	✓
Safety Valve	✓	✓	✓	✓
Low Water Cut Off	✓	✓	✓	✓
Moisture Trap	✓	✓	✓	✓
Dial Thermometer	✓	✓	✓	✓
Analog Pressure Gauge	✓	✓	✓	✓
Analog Compound Gauge	✓	✓	✓	✓
Automatic Pressure Control	✓	✓	✓	✓
Automation Device	●	DTC	PID	PLC
Automatic Ster Controlled Valves	●	✓	✓	✓
Digital Temperature Controller	●	✓	✓	✓
Digital Temperature Display	●	✓	✓	✓
Digital Timer	●	✓	✓	✓
Alarm	●	✓	✓	✓
Digital Pressure Display	●	●	✓	✓
Printer	●	●	optional	✓
Touch Screen Display	●	●	●	✓
Pre Set Sterilization Cycle	●	●	●	✓
Automatic Water feeding	●	●	●	✓
Auto Exhaust	●	●	●	✓

Customization of products can be done as per client's requirements

Our Esteemed Clients

Hospitals and Medical Colleges



Pharmaceutical Lab



Food and Beverage



Agricultural Labs



Get in Touch

Marketed by



📍 : D-149 Preet Vihar
New Delhi 110092.
✉ : ysupltd@gmail.com
🌐 : yorcoscientific.com
☎ : 8595293382

Manufactured by



YORCO Scientific Udyog Pvt.Ltd
📍 : 13 Loni road, site 2,
Mohan nagar, Ghaziabad,
Uttar Pradesh 201007.
✉ : ysupltd@gmail.com
🌐 : yorcoscientific.com
☎ : 935000607